



Roho® Low Profile® Quadro® Select Cushion

Manufacturer:
The Roho

The Low Profile Quadro Select Cushion sets the standard in the industry for overall performance in wheelchair seating. Each and every user will benefit from a customized fit through the simple push of a knob. The revolutionary Isoflo® Memory Control offers shape shifting capabilities while the user is seated, allowing quick and easy on demand adjustment to maximize function. With the built-in stability and simplicity of the Low Profile Quadro Select Cushion, no longer will you have to sacrifice maximum skin protection to get stability, positioning or convenience. The Low Profile Quadro Select Cushion does it all.

Roho® Low Profile® Quadro® Select Cushion

FEATURES:

The Roho



- **Positioning:** The Quadro Select Cushion allows for quick, simple position selection for every user while seated. It allows the user to selectively change their position throughout the day as their particular activity demands.
- **Adjustable:** Fits specific seating needs and unique body shapes.
- **Skin Protection:** The Roho Dry Floatation® technology provides unparalleled performance in facilitating blood flow by fitting, matching and tracking the shape of a person.
- **Stability:** The Low Profile Quadro Select Cushion offers a new level of stability never before found in an air flotation cushion.

SPECIFICATIONS

LIGHTWEIGHT	Easy to handle and carry, at just 2.5 lbs/1.1 kg (QS99LPC). Varies by size.
WASHABLE	Cushion can be easily cleaned with soap and water. Cover is machine washable.
REPAIRABLE	Neoprene rubber cushion can be patched to extend use.
VENTILATED	Separated air cells allow air circulation to reach the surface of the skin and provide a cooling environment.
INCLUDES	Low Profile Quadro Select Cushion, two-way stretch cover, hand inflation pump, repair kit, and operations manual.

CONSTRUCTION	Flame Resistant. Neoprene Rubber.
SIZES	Wide range to fit chair widths from 12" to 20"/30.5 cm to 51cm
WARRANTY	3 Year Limited Warranty