

With features like TB-Flex seating, Mid-Wheel 6® Drive Design and ATX Suspension, the Q4 is engineered for active users. The Q4 delivers exceptional rehab capability and performance at an excellent value.



Specifications

DRIVE WHEELS	Standard: 14" Solid Optional: 14" Pneumatic
CASTER WHEELS	Front: 6" Solid Rear: 6" Solid
ANTI-TIP WHEELS	N/A
MAXIMUM SPEED ¹	Up to 6 mph
GROUND CLEARANCE	2.5" (at caster pivot)
TURNING RADIUS ²	20"
OVERALL LENGTH ²	35.5"
BASE WIDTH	24"
SEATING SIZES	
TB-Flex Seat	Width: 15-22" Depth: 15-20"
TRU-Balance® 2 Power Tilt	Width: 16-20" Depth: 16-20"
SEAT-TO-FLOOR HEIGHTS	
TB-Flex Seat ⁵	15.875" - 17.875"
TRU-Balance® 2 Power Tilt	17.3" - 19.3"
MANUAL RECLINE 6	Yes
BATTERY SIZE ³	NF-22
BATTERY WEIGHT ⁷	38 lbs.
AVAILABLE ELECTRONICS 8	120A Q-Logic 2 EX 75A Q-Logic 2 NE+
BATTERY CHARGER	8A Off-board
MOTOR PACKAGES	4-Pole Mid-Wheel 6®
WEIGHT CAPACITY ⁴	300 lbs.
BASE WEIGHT	138.6 lbs.
	*

Features

- Standard 6 mph motors offer advanced performance, precise control and responsiveness
- Bold aesthetics
- Q-Logic 2 NE+ electronics standard
- High-back seat
- ATX Suspension (Active-Trac® with extra stability) incorporates front OMNI-Casters and semi-independent rear caster beam for enhanced performance over more varied terrain

Options

- Q-Logic 2 EX Drive Controls
- Accu-Trac Advanced Tracking Technology
- TB-Flex seating system
- TRU-Balance® 2 Power Tilt
- High-back seat with solid seat pan

Colors







Red Corvette

Electric Blue

in Black

- 1) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.
- 2) Without front riggings
- 3) Battery capacity may change when adding power positioning systems.
- 4) Weight capacity is dictated by the selection of seating.
 5) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.
- 6) Maximum dimensions of the Manual Recline vary by power base and seating system. Call for details.
- 7) Battery weight may vary +/- 3 lbs. based on manufacturer.
- 8) All Quantum Rehab® units can be ordered with specialty controls.